

# Prestained Protein Ladder, (10-170kdA)

Cat.No. SL7011(PR911654) Size:500µl

Cat.No. SL7012) Size:100µl

- 3µl or 5µl per loading for clear visualization during electrophoresis on 15-well or 10-well mini-gel, respectively.
- 2 ~ 3 μl per well for general Western transferring.
- Apply more for thicker (> 1.5 mm) or larger gel.

## **Descripition:**

The prestained protein ladder is a three-color protein standard with 10 pre-stained proteins covering a wide range molecular weights for 10 to 180 KDa. Proteins are covalently coupled with a blue chromophore except for two reference bands(one green and one red band at 25 kDa and 75 kDa respectively) when separated on SDS-PAGE (Tris-Glycine buffer) .The Prestained protein ladder is designed for monitoring protein separated during SDS polyacrylamide gel electrophoresis , Verification of Western transfer efficiency on membranes (PVDF , nylon or nitrocellulose) and for approximate sizing of proteins. The Ladder is Supplied in gel loading buffer and is ready to use. Do not heat, dilute , add reducing agent before loading.

### **Contents:**

Approximately 0.2  $\sim$  0.4 mg/ml of each protein in buffer (20 mM Trisphosphate pH 7.5 at 25  $^{\circ}$ C), 2% SDS, 1 mM 2-Mercaptoethanol, 3.6 M Urea, and 15% (v/v) Glycerol).

## **Quality Control:**

 $5~\mu$  of Prestained protein Ladder resolves 10 bands in 4-20% SDS-PAGE (Tris-Glycine buffer) and after Western blotting to PVDF membrane.

### Storage:

Stable for up to 2 weeks at 25 °C.

Note. The apparent molecular weight of each protein (KDa) has been determined by calibration against an unstained protein ladder in each electrophoresis condition.

\*supplement data should be considered for more accurate adjustment.

All product are for research use only. Caution: Not intended for human or animal diagnostic or therapeutic uses.

#### SinaClon BioScience


www.sinaclon.com

Tel: +98(0)21 4463 0050-51 +98(0)21 4466 5156

Customer Service: +98(0)21 4463 3016 Place your order: <u>order@sinaclon.com</u> Guide for Molecular Weight Estimation(KDa):
Migration patterns of Prestained Protein Ladder in
different electrophoresis conditions are listed below.

Stable for up to 3 months at 4 °C. For long term storage

at -20 °C.



Sale Office: No. 16, Sayeh Building, Nafisi Ave., Phase 1, Ekbatan, Tehran, Iran

Tell: +9821 44661950, +9821 44661903